

VITAPAN PLUS – A symbiosis of function and aesthetics

On the basis of the popular classic VITAPAN, VITA Zahnfabrik has now developed the compact anterior tooth assortment VITAPAN PLUS. Taking into account, and strictly adhering to the wishes of customers from all over the world, the result is a modern anterior tooth line, which thanks to convincing interior and exterior characteristics, is suitable for use by every prosthodontist, and for every indication. VITAPAN PLUS shall be available as from May 2011, matched to the VITA SYSTEM 3D-MASTER and VITAPAN Classical A1-D4.

The new VITAPAN PLUS forms the perfect symbiosis of function and aesthetics. In this way, thanks to the authentic surface morphology and mamelon structure and balanced opalescence, luminescence and translucency, the light reflection properties exhibit an extremely lifelike and richly nuanced play of colours. Details such as the whitish incisal margins, pronounced angular characteristics and a type-appropriate labial marginal ridge design likewise contribute to the overall natural aesthetics. The horizontal and vertical curvature characteristics were featured in the design with a view to effective lip support to conserve the patient's individual physiognomy and to enable good phonetics, and the nature-identical design of the palatal surfaces ensures the accustomed tongue sensation. The slightly broadened tooth necks with a stepless transition to the tooth crown furthermore allow patient-oriented and hygiene-friendly gingival modelling with flexible design of the interdental spaces.

The innovative tooth line is suitable both for complete denture and implant prosthetics as well as attachments, telescopic and conus



constructions and model casting prosthesis. VITAPAN PLUS is reliable and accurate to process and demonstrates extraordinary aesthetics and reliable wear characteristics. The homogeneous material with excellent strength values guarantees a long life of the restorations.

Furthermore, the VITAPAN Plus anterior tooth can be combined with all VITA posterior tooth lines.

For further information visit our website www.vita-zahnfabrik.com or call the VITA hotline on +49 (0)7761-562 222.



Fig. 1: VITAPAN PLUS – The modern anterior tooth line for every prosthodontist and every indication.

We would appreciate a sample copy of any reprints.

VITA Zahnfabrik H. Rauter GmbH & Co. KG

Headquartered in Bad Säckingen/Germany, VITA Zahnfabrik H. Rauter GmbH & Co. KG has been developing, producing and marketing innovative solutions for dental prosthetics according to consistently high quality standards for over 85 years, and has been known from the very beginning as a pioneer and worldwide trendsetter. The VITA shade standard, for instance, is recognised internationally in the dental branch as a shade reference system. Users in 120 different countries benefit from the comprehensive range of products and services provided by VITA Zahnfabrik. These include analogue and digital tooth shade determination systems, acrylic and ceramic teeth, veneering and framework materials for



VITA

conventional and computer-aided manufacturing procedures, dental equipment as well as a wide range of service and training facilities.

Your contact person:

Gert Schimak

Tel. + 49 (0) 7761 562 226

Fax + 49 (0) 7761 562 281

E-Mail: g.schimak@vita-zahnfabrik.com

Internet: www.vita-zahnfabrik.com

